

Series	Gas Type	Product Name	Life	Maximum Overload Range	Details	Fixed/Portable	Typical Applications
CiTiceL [®]	O ₂	C/NLH	Proportional	0-1000ppm	Metal can sensor for trace O ₂ with low H ₂ cross sensitivity	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	C/NLL	Proportional	0-1000ppm	Metal can sensor for trace O ₂ - high output	Fixed/Portable	Industrial safety
5 Series	O ₂	5FO	2 years in air	0-30%	Flue gas sensor	Fixed/Portable	Emissions
7 Series	O ₂	7OX-V	2 years in air	0-30%	Compact sensor - plastic can. For use in extreme environments	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	2FO	2 years in air	0-30%	Plastic can sensor	Fixed/Portable	Emissions
CiTiceL [®]	O ₂	2FO-N	2 years in air	0-30%	Plastic can sensor with capillary nose	Fixed/Portable	Emissions
CiTiceL [®]	O ₂	C/2	18 months in air	0-30%	Half size metal can sensor	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	C/2PN	1 yr in air	0-30%	Metal can sensor with capillary nose	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	C/N	9 months in air	0-30%	Metal can sensor with capillary nose	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	C/S	9 months in air	0-30%	Metal can sensor	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	C/Y	1 year in air	0-30%	General purpose metal can sensor	Fixed/Portable	Industrial safety
7 Series	O ₂	T7OX-V	2 years in air	0-25%	Compact sensor - plastic can with transmitter board	Fixed/Portable	Industrial safety
CiTiceL [®]	O ₂	4OXV	2 years in air	0-25%		Fixed/Portable	Industrial safety
AutoO ₂ [®]	O ₂	AO2	2 years in air	0-100%	Automotive sensor with molex connector - BAR97/PTB Approved	Fixed/Portable	Automotive
AutoO ₂ [®]	O ₂	AO3	2 years in air	0-100%	Automotive sensor with 0.141" miniature phone jack socket - BAR97/PTB Approved	Fixed/Portable	Automotive
MediceL [®]	O ₂	MOX-1	20 Months in 100% O ₂	0-100%	Sensor with Molex connector, M16 thread	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-2	20 Months in 100% O ₂	0-100%	Sensor with jack socket connector, M16 thread	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-3	20 Months in 100% O ₂	0-100%	Sensor with UK telephone socket, M16 thread	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-4	20 Months in 100% O ₂	0-100%	Teledyne T-7 compatible sensor	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-6	14 Months in 100% O ₂	0-100%	Single cathode Draeger compatible sensor, slip rings	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-9	14 Months in 100% O ₂	0-100%	Sensor with Molex Connector, M16 Thread, Siemens compatible	Fixed/Portable	Medical - Anaesthesia agent compatible
MediceL [®]	O ₂	MOX-20	2 years in air	0-100%	Fast response sensor, Molex connector, M16 thread	Fixed/Portable	Medical - Breath by breath

MediceL [®]	O ₂	Divecel 3	2 years in air	0-100%	Standard output sensor, Molex connector, M16 thread	Fixed/Portable	Medical - Diving
MediceL [®]	O ₂	DO2	2 years in air	0-100%	Double output sensor, Molex connector, M16 thread	Fixed/Portable	Medical - Diving
MediceL [®]	O ₂	InQOX	2 years in 50% O ₂	0-100%	Neonatal incubator sensor, Molex connector, M16 thread	Fixed/Portable	Medical - Incubators